

Anti-Geminin Rabbit antibody



Catalog Number: 168960

Product name

Anti-Geminin Rabbit antibody

Specificity

Human

Antibody description

Geminin Rabbit monoclonal antibody

Preparation

Antigen: Recombinant protein of human Geminin

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB: 1/1000

Validations



Western blot detection of Geminin in Jurkat cell lysates using Geminin Rabbit mAb(1:1000 diluted).Predicted band size:24kDa.Observed band size:24kDa.