

H6N2 HA (His Tag) recombinant protein



Catalog Number: 501260

General Information

Protein Construction

A DNA sequence encoding the Influenza A virus (A/chicken/Guangdong/C273/2011(H6N2)) hemagglutinin (AFJ19048.1) (Met1-Arg344), termed as HA1, was expressed with a polyhistidine tag at the C-terminus.

Organism

H6N2

Expression Host

Human Cells

QC Testing

Purity

> 95 % as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Asp 17

Molecular Mass

The recombinant HA1 subunit of the Influenza A virus (A/chicken/Guangdong/C273/2011(H6N2))

hemagglutinin consists 339 amino acids and predicts a molecular mass of 38.4 kDa.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Influenza A H6N2 (A/chicken/Guangdong/C273/2011) Hemagglutinin / HA1 Protein (His Tag)