

# SH3GL3 antibody



Catalog Number: 115260

## Product name

SH3GL3 antibody

## Immunogen

[Human SH3GL3 Recombinant protein \(GST tag\)](#)

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

SH3GL3 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue. Observed molecular weight by Western-blot: 39.3kd

## Preparation

This antibody was obtained by immunization of SH3GL3 recombinant protein (Accession Number: BC042864). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:200-1:2000

## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:115260(SH3GL3 antibody) at dilution of 1:800