

# ZO-1 antibody

Catalog Number: 117260

---

## Product name

ZO-1 antibody

## Specificity

Human; other species not tested.

## Antibody description

ZO-1 Rabbit Polyclonal antibody. Positive IF detected in HUVEC cells.

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_003257).  
Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol  
pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IF

## Dilutions

Recommended Dilution:

IF: 1:20-1:200

## Validations



Immunofluorescent analysis of HUVEC cells using Catalog No:117260(TJP1 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)