

PDE3B antibody



Catalog Number: 113660

Product name

Polyclonal

PDE3B antibody

Ig Type

Specificity

Rabbit IgG

Human; other species not tested.

Applications

Antibody description

ELISA, IP

PDE3B Rabbit Polyclonal antibody. Positive IP detected in HeLa cells. Observed molecular weight by Western-blot: 124 kDa

Dilutions

Recommended Dilution:

Preparation


IP: 1:200-1:1000

This antibody was obtained by immunization of PDE3B recombinant protein (Accession Number: NM_000922). Purification method: Antigen affinity purified.

Validations

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

IP Result of anti-PDE3B (IP:Catalog No:113660, 5ug; Detection:Catalog No:113660 1:300) with HeLa cells lysate 800ug.

Clonality