

# Anti-Caspase-3 Rabbit antibody



Catalog Number: 166660

## Product name

Anti-Caspase-3 Rabbit antibody

## Specificity

Human, Mouse

## Antibody description

Caspase-3 Rabbit polyclonal antibody

## Preparation

Antigen: Purified recombinant human Caspase-3 protein fragments expressed in E.coli.

## Formulation

PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB


## Dilutions

WB: 1/1000

## Validations



Western blot detection of Caspase-3 in Jurkat and 3T3 cell lysates using Caspase-3 Rabbit pAb (1:1000 diluted). Predicted band size: 35KDa.



Observed band size:35KDa.

Western blot analysis of extracts from Jurkat cells, untreated (-) or treated with Staurosporine (1  $\mu$ M; 3 hr), using Caspase-3 Rabbit pAb (166660, upper) or  $\beta$ -Actin (8F10) Mouse mAb (lower).