

# Anti-Sox10 antibody

Catalog Number: 176660

## Product name

Anti-Sox10 antibody

## Specificity

Mouse

## Antibody description

Mouse monoclonal antibody to Sox10

## Preparation

This antigen of this antibody was recombinant protein

## Formulation

Liquid, 1\*PBS (pH7.4), 0.2% BSA, 40% Glycerol.  
Preservative: 0.05% Sodium Azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications

WB, IHC-P

## Dilutions

WB: 1:500-1:1000

IHC-P: 1:50

## Validations





Fig1: Western blot analysis on F9 (1) and mouse fetal brain (2) cell lysates using anti-SOX10 Mouse



mAb.

Fig2: Immunohistochemical analysis of paraffin-embedded mouse fetal brain (13.5 days) tissue using anti- SOX10 Mouse mAb.