

# Anti-FNTB Rabbit antibody



Catalog Number: 168660

## Product name

Rabbit IgG

Anti-FNTB Rabbit antibody

## Applications

## Specificity

WB, FC, IP

Human

## Dilutions

## Antibody description

WB: 1/2000

FNTB Rabbit monoclonal antibody

FC: 1/20

## Preparation

IP: 1/20

Antigen: A synthetic peptide of human FNTB

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage


Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

## Validations



Western blot detection of FNTB in K562, HeLa cell lysates using FNTB Rabbit mAb(1:1000 diluted). Predicted band size: 49kDa. Observed band size: 49kDa.