

Human UBE2R2/CDC34B/UBC3B (His tag) recombinant protein



Catalog Number: 519660

General Information

Gene Name Synonym

Ubiquitin-Conjugating Enzyme E2 R2, Ubiquitin Carrier Protein R2, Ubiquitin-Conjugating Enzyme E2-CDC34B, Ubiquitin-Protein Ligase R2, UBE2R2, CDC34B, UBC3B

Protein Construction

A DNA sequence encoding the Human UBE2R2/CDC34B/UBC3B Met1-Ser238 is expressed with a 6His tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 50mm HEPES, 150mM NaCl, pH 7.5.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human UBE2R2/CDC34B/UBC3B (His tag)
recombinant protein