

P2RX3 Polyclonal Antibody



Catalog Number: 160360

Product name

P2RX3 Polyclonal Antibody

Specificity

Rat

Antibody description

Polyclonal antibody to P2RX3

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 44-317 of human P2RX3 (NP_002550.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:1000 - 1:2000

Validations



Western blot - P2RX3 Polyclonal Antibody

Western blot analysis of extracts of rat testis, using P2RX3 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 20s.