

Human ALK-3 / BMPR1A (His & Fc Tag) recombinant protein



Catalog Number: 500360

General Information

Gene Name Synonym

Activin receptor-like kinase 3; Serine/threonine-protein kinase receptor R5

Protein Construction

A DNA sequence encoding the extracellular domain (Met 1-Arg 152) of human ALK3 (NP_004320.2) was fused with the C-terminal polyhistidine-tagged Fc region of human IgG1 at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Activity

Measured by its ability to inhibit BMP4 induced activity in MC3T3-E1 Mouse osteoblastic cells. The ED₅₀ for this effect is typically 0.1-0.3 µg/ml in the presence of 50 ng/mL of recombinant human BMP4.

Purity

> 97 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per µg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Gln 24

Molecular Mass

The recombinant human ALK3/Fc is a disulfide-linked homodimeric protein after removal of the signal peptide. The reduced monomer consists of 376 amino acids and has a predicted molecular mass of 42 kDa. In SDS-PAGE under reducing conditions, the apparent molecular mass of rh ALK3/Fc monomer is approximately 56 kDa due to glycosylation.

Formulation

Lyophilized from sterile PBS, pH 7.4
1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human ALK-3 / BMPR1A (His & Fc Tag) recombinant protein



Catalog Number: 500360

Human BMPRIA / ALK-3 / CD292 Protein (His & Fc Tag) SDS-PAGE