

# HIST1H2AC antibody



Catalog Number: 111360

## Product name

HIST1H2AC antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

HIST1H2AC Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue, rat liver tissue. Observed molecular weight by Western-blot: 17 kDa

## Preparation

This antibody was obtained by immunization of HIST1H2AC recombinant protein (Accession Number: BC017379). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

## Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:111360(HIST1H2AC Antibody) at dilution of 1:300