

# Anti-Tnni3 antibody



Catalog Number: 176360

## Product name

Anti-Tnni3 antibody

## Specificity

Mouse, Rat

## Antibody description

Rabbit polyclonal antibody to Tnni3

## Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from mouse tnni3. aa range:5-54

## Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, IF

## Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

## Validations




Fig1:: Western Blot analysis of mouse cells using Troponin I-C Polyclonal Antibody. Secondary antibody (catalog#: HA1001) was diluted at 1:20000