

KCNN4 antibody



Catalog Number: 107360

Product name

KCNN4 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

KCNN4 Mouse Monoclonal antibody. Positive WB detected in human blood tissue. Observed molecular weight by Western-blot: 48 kDa

Preparation

This antibody was obtained by immunization of KCNN4 recombinant protein (Accession Number: NM_002250). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2a

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human blood tissue were subjected to SDS PAGE followed by western blot with Catalog No:107360(KCNN4 Antibody) at dilution of 1:1000