

Human FANCA / FACA (His Tag) recombinant protein



Catalog Number: 502760

General Information

Protein Construction

A DNA sequence encoding the human FANCA (AAH08979.1) (Met 1-Cys 297) was expressed, with a polyhistidine tag at the N-terminus.

Organism

Human

Expression Host

Baculovirus-Insect Cells

QC Testing

Purity

> 85 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Met

Molecular Mass

The recombinant human FANCA consists of 316 amino acids and predicts a molecular mass of 35.2 kDa. It migrates as an approximately 35 kDa band

in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile 20mM Tris, 200mM NaCl, pH 8.5, 20% gly, 0.1% Triton, 1mM TCEP

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide

Storage


Store it under sterile conditions at -20°C to -80°C .

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human FANCA / FACA Protein (His Tag) SDS-PAGE