

SYK antibody

Catalog Number: 115760

Product name

SYK antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

SYK Rabbit Polyclonal antibody. Positive IHC detected in human spleen.

Preparation

This antibody was obtained by immunization of SYK recombinant protein (Accession Number: XM_005252147). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, IHC


Dilutions

Recommended Dilution:

IHC: 1:20-1:200

Validations

Immunohistochemical of paraffin-embedded human spleen using Catalog No:115760(SYK antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human spleen using Catalog No:115760(SYK antibody) at dilution of 1:50 (under 40x lens)