

CACNA1S antibody



Catalog Number: 108760

Product name

CACNA1S antibody

Specificity

Human; other species not tested.

Antibody description

CACNA1S Rabbit Polyclonal antibody. Positive WB detected in human brain tissue. Positive IHC detected in human skeletal muscle tissue. Observed molecular weight by Western-blot: 200-220kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_000069). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:

WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:108760(CACNA1S antibody) at dilution of



1:1000
Immunohistochemistry of paraffin-embedded human skeletal muscle slide using Catalog No:108760(CACNA1S Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human skeletal muscle slide using Catalog No:108760(CACNA1S Antibody) at dilution of 1:50 (under 40x lens)