

SARDH antibody



Catalog Number: 115070

Product name

SARDH antibody

Specificity

Human, Rat; other species not tested.

Antibody description

SARDH Rabbit Polyclonal antibody. Positive WB detected in rat liver tissue. Observed molecular weight by Western-blot: 100 kDa

Preparation

This antibody was obtained by immunization of SARDH recombinant protein (Accession Number: NM_001134707). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



rat liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:115070(SARDH Antibody) at dilution of 1:1500