

Anti-MED8 antibody



Catalog Number: 176070

Product name

Anti-MED8 antibody

Specificity

Human, Mouse, Rat

Antibody description

Rabbit polyclonal antibody to MED8

Preparation

This antigen of this antibody was synthesized peptide derived from med8 . at aa range: 70-150

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations




Fig1.; Western Blot analysis of JK cells using Med8 Polyclonal Antibody diluted at 1: 500 cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Inventbiotech, MN, USA).