

Human LCN2/HNL/NGAL (His tag) recombinant protein



Catalog Number: 518070

General Information

Gene Name Synonym

Neutrophil gelatinase-associated lipocalin, NGAL, 25 kDa alpha-2-microglobulin-related subunit of MMP-9, Lipocalin-2, Oncogene 24p3, Siderocalin LCN2, p25, HNL, NGAL

Protein Construction

A DNA sequence encoding the Human LCN2/HNL/NGAL Gln21-Gly198 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of PBS, 50% glycerol, pH7.4.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human LCN2/HNL/NGAL (His tag) recombinant protein