

# MAOB antibody



Catalog Number: 112470

## Product name

MAOB antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

MAOB Rabbit Polyclonal antibody. Positive IHC detected in human gliomas tissue. Positive WB detected in human liver tissue, mouse liver tissue. Observed molecular weight by Western-blot: 59 kDa

## Preparation

This antibody was obtained by immunization of MAOB recombinant protein (Accession Number: NM\_000898). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

## Validations



human liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:112470(MAOB antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:112470(MAOB antibody) at dilution of 1:100 (under 10x lens)