

BMSC UbP antibody



Catalog Number: 108470

Product name

BMSC UbP antibody

WB: 1:500-1:5000

IP: 1:500-1:5000

Specificity

Human, Mouse, Rat; other species not tested.

IHC: 1:20-1:200

IF: 1:10-1:100

Antibody description

BMSC UbP Rabbit Polyclonal antibody. Positive IHC detected in human gliomas tissue. Positive IF detected in HepG2 cells. Positive IP detected in HeLa cells. Positive WB detected in HepG2 cells, HeLa cells, Jurkat cells, mouse testis tissue, rat testis tissue. Observed molecular weight by Western-blot: 41 kDa

Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:117205(UBL7 antibody) at dilution of 1:400

Preparation

This antibody was obtained by immunization of BMSC UbP recombinant protein (Accession Number: NM_032907). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications


ELISA, WB, IHC, IP, IF

Dilutions

Recommended Dilution:




IP Result of anti-UBL7 (IP:Catalog No:117205, 4ug; Detection:Catalog No:117205 1:1000) with HeLa cells lysate 2000ug.



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:117205(UBL7 antibody) at dilution of 1:50 (under 10x lens)

BMSC UbP antibody

Catalog Number: 108470



Immunofluorescent analysis of HepG2 cells using Catalog No:117205(UBL7 Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)