

## General Information

### Protein Construction

A DNA sequence encoding the Influenza A virus (A/turkey/Germany-MV/R2472/2014 (H5N8)) hemagglutinin (translated amino acids of EPI544756) (Met1-Gln530), termed as HA, was expressed with a polyhistidine tag at the C-terminus.

### Organism

H5N8

### Expression Host

Baculovirus-Insect Cells

## QC Testing

### Purity

> 90 % as determined by SDS-PAGE.

### Endotoxin

< 1.0 EU per  $\mu\text{g}$  protein as determined by the LAL method.

### Stability

Samples are stable for up to twelve months from date of receipt at  $-70^{\circ}\text{C}$

### Predicted N terminal

Asp 17

### Molecular Mass

The recombinant hemagglutinin of Influenza A virus (A/turkey/Germany-MV/R2472/2014 (H5N8))

consists of 525 amino acids and predicts a molecular mass of 60.1 kDa.

### Formulation

Lyophilized from sterile 20 mM Tris, pH 8.0, 500 mM NaCl, 10 % gly.

1. Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization. Specific concentrations are included in the hardcopy of COA.
2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$ . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

A hardcopy of COA with reconstitution instruction is sent along with the products. Please refer to it for detailed information.

### SDS-PAGE



H5N8 HA Protein 14130