

TBCB antibody



Catalog Number: 115870

Product name

TBCB antibody

IP: 1:200-1:2000

IF: 1:10-1:100

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TBCB Rabbit Polyclonal antibody. Positive IP detected in HeLa cells. Positive WB detected in HeLa cells, human brain tissue. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 32 kDa

Preparation

This antibody was obtained by immunization of TBCB recombinant protein (Accession Number: NM_001281). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IF, IP

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:115870(TBCB antibody) at dilution of 1:500



Immunofluorescent analysis of HeLa cells, using TBCB antibody Catalog No:115870 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).



IP Result of anti-TBCB (IP:Catalog No:115870, 3ug; Detection:Catalog No:115870 1:500) with

TBCB antibody



Catalog Number: 115870

HeLa cells lysate 1200ug.