

# TBCB antibody

Catalog Number: 115870

## Product name

TBCB antibody

IP: 1:200-1:2000

IF: 1:10-1:100

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

TBCB Rabbit Polyclonal antibody. Positive IP detected in HeLa cells. Positive WB detected in HeLa cells, human brain tissue. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 32 kDa

## Preparation

This antibody was obtained by immunization of TBCB recombinant protein (Accession Number: NM\_001281). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB, IF, IP

## Dilutions


Recommended Dilution:

WB: 1:500-1:5000


## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:115870(TBCB antibody) at dilution of 1:500



Immunofluorescent analysis of HeLa cells, using TBCB antibody Catalog No:115870 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).



IP Result of anti-TBCB (IP:Catalog No:115870, 3ug; Detection:Catalog No:115870 1:500) with

# TBCB antibody



Catalog Number: 115870

---

HeLa cells lysate 1200ug.