

TBCB antibody

Catalog Number: 115870

Product name

TBCB antibody

IP: 1:200-1:2000

IF: 1:10-1:100


Specificity

Human, Mouse, Rat; other species not tested.

Validations

Antibody description

TBCB Rabbit Polyclonal antibody. Positive IP detected in HeLa cells. Positive WB detected in HeLa cells, human brain tissue. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 32 kDa




Preparation

This antibody was obtained by immunization of TBCB recombinant protein (Accession Number: NM_001281). Purification method: Antigen affinity purified.

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:115870(TBCB antibody) at dilution of 1:500

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Immunofluorescent analysis of HeLa cells, using TBCB antibody Catalog No:115870 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).

Ig Type

Rabbit IgG


Applications

ELISA, WB, IF, IP

Dilutions

Recommended Dilution:

WB: 1:500-1:5000



IP Result of anti-TBCB (IP:Catalog No:115870, 3ug; Detection:Catalog No:115870 1:500) with

TBCB antibody



Catalog Number: 115870

HeLa cells lysate 1200ug.