

# APOB Monoclonal Antibody



Catalog Number: 166570

## Product name

APOB Monoclonal Antibody

## Specificity

Human

## Antibody description

Monoclonal antibody to APOB

## Preparation

Antigen: Recombinant protein of human APOB

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Monoclonal

## Ig Type

Mouse IgG

## Applications

WB, ICC, FC


## Dilutions

WB 1:500 - 1:2000

ICC 1:500 - 1:1000

FC 1:50 - 1:200

## Validations



## Western blot - APOB Monoclonal Antibody

Western blot analysis of extracts of Recombinant protein, using APOB antibody. Secondary antibody: HRP Goat Anti-Mouse IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.