

Anti-NOTCH4 Rabbit antibody



Catalog Number: 167570

Product name

Anti-NOTCH4 Rabbit antibody

Specificity

Human

Antibody description

NOTCH4 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human NOTCH4

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB: 1/2000-1/10000

Validations



Western blot detection of NOTCH4 in HeLa,U87-MG cell lysates using NOTCH4 Rabbit mAb(1:1000 diluted).Predicted band size:36kDa.Observed band size:36kDa.