

Human NT5E/NT5/NTE (His tag) recombinant protein



Catalog Number: 518570

General Information

Gene Name Synonym

5'-Nucleotidase, 5'-NT, Ecto-5'-Nucleotidase, CD73, NT5E, NT5, NTE

Protein Construction

A DNA sequence encoding the Human NT5E/NT5/NTE Trp27-Lys547 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM TrisHCl, 120mM NaCl, 4mM CaCl₂, 20% Glycerol, pH 7.5.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human NT5E/NT5/NTE (His tag) recombinant protein