

# DNMT3L antibody



Catalog Number: 109970

---

## Product name

DNMT3L antibody

## Specificity

human; other species not tested.

## Antibody description

DNMT3L Rabbit Polyclonal antibody. Positive WB detected in PC-3 cells, human testis tissue. Observed molecular weight by Western-blot: 44kd and 70kd

## Preparation

This antibody was obtained by immunization of DNMT3L recombinant protein (Accession Number: NM\_175867). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:109970(DNMT3L antibody) at dilution of 1:300