

HBM-Specific antibody



Catalog Number: 111270

Product name

HBM-Specific antibody

Specificity

Human; other species not tested.

Antibody description

HBM-Specific Rabbit Polyclonal antibody. Positive WB detected in human heart tissue, recombinant protein. Observed molecular weight by Western-blot: 16 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_005332). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions


Recommended Dilution:

WB: 1:200-1:2000

Validations



WB results of Catalog No:111270 (HBM antibody)



at dilution of 1:2000.

human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:111270(HBM-Specific antibody) at dilution of 1:300