

NR5A2 antibody



Catalog Number: 113270

Product name

NR5A2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

NR5A2 Rabbit Polyclonal antibody. Positive WB detected in mouse ovary tissue, rat testis tissue. Observed molecular weight by Western-blot: 61-67 kDa

Preparation

This antibody was obtained by immunization of NR5A2 recombinant protein (Accession Number: NM_003822). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



mouse ovary tissue were subjected to SDS PAGE followed by western blot with Catalog No:113270(NR5A2 Antibody) at dilution of 1:1000