

Anti-PPP4R1L antibody



Catalog Number: 176270

Product name

Anti-PPP4R1L antibody

Specificity

Human, Monkey

Antibody description

Rabbit polyclonal antibody to PPP4R1L

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human ppp4r1l. aa range:191-240

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations




Fig1.; Western Blot analysis of various cells using PPP4R1L Polyclonal Antibody