

FKBP14 antibody



Catalog Number: 110670

Product name

FKBP14 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

FKBP14 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, HeLa cells, human brain tissue. Observed molecular weight by Western-blot: 28 kDa

Preparation

This antibody was obtained by immunization of FKBP14 recombinant protein (Accession Number: NM_017946). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:110670(FKBP14 antibody) at dilution of 1:500