

PDE8A antibody

Catalog Number: 113670

Product name

PDE8A antibody

IP: 1:500-1:5000

IHC: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PDE8A Rabbit Polyclonal antibody. Positive IP detected in mouse testis tissue. Positive WB detected in mouse testis tissue. Positive IHC detected in human colon cancer tissue. Observed molecular weight by Western-blot: 93kd, 88kd

Preparation

This antibody was obtained by immunization of PDE8A recombinant protein (Accession Number: NM_002605). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IP

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:113670(PDE8A antibody) at dilution of 1:1000



IP Result of anti-PDE8A (IP:Catalog No:113670, 5ug; Detection:Catalog No:113670 1:1000) with mouse testis tissue lysate 8000ug.




Immunohistochemical of paraffin-embedded human colon cancer using Catalog

PDE8A antibody

Catalog Number: 113670

No:113670(PDE8A antibody) at dilution of 1:100
(under 10x lens)



Immunohistochemical of paraffin-embedded
human colon cancer using Catalog
No:113670(PDE8A antibody) at dilution of 1:100
(under 40x lens)