

# UBXN11 antibody



Catalog Number: 116670

## Product name

UBXN11 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

UBXN11 Rabbit Polyclonal antibody. Positive IHC detected in human lung cancer tissue, human testis tissue. Positive IP detected in mouse testis tissue. Positive WB detected in PC-3 cells, mouse testis tissue. Observed molecular weight by Western-blot: 48-54kd

## Preparation

This antibody was obtained by immunization of UBXN11 recombinant protein (Accession Number: BC038106). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC, IP

## Dilutions


Recommended Dilution:

WB: 1:200-1:2000


IP: 1:200-1:2000

IHC: 1:20-1:200

## Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:116670(UBXN11 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human lung cancer using Catalog No:116670(UBXN11 antibody) at dilution of 1:100 (under 10x lens)

IP Result of anti-UBXN11 (IP:Catalog No:116670, 4ug; Detection:Catalog No:116670 1:500) with mouse testis tissue lysate 4000ug.