

Anti-Caspase-11 Rabbit antibody



Catalog Number: 168670

Product name

Anti-Caspase-11 Rabbit antibody

Ig Type

Rabbit IgG

Specificity

Mouse

Applications

WB, IP

Antibody description

Caspase-11 Rabbit monoclonal antibody

Dilutions

1:1000

Preparation


Antigen: Recombinant protein of mouse Caspase-11

1:20

Validations

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.



Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Western blot detection of Caspase-11 in C6,3T3,Hela cell lysates using Caspase-11 Rabbit mAb(1:1000 diluted).Predicted band size:43kDa.Observed band size: 43,38,25kDa.

Clonality

Monoclonal