

C15orf53 antibody



Catalog Number: 108670

Product name

C15orf53 antibody

Specificity

Human; other species not tested.

Antibody description

C15orf53 Rabbit Polyclonal antibody. Positive WB detected in A431 cells, COLO 320 cells, HepG2 cells. Positive IHC detected in human tonsillitis tissue, human spleen tissue. Observed molecular weight by Western-blot: ~30 kDa

Preparation

This antibody was obtained by immunization of C15orf53 recombinant protein (Accession Number: NM_207444). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



Immunohistochemistry of paraffin-embedded human tonsillitis slide using Catalog No:108670(C15orf53 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human tonsillitis slide using Catalog No:108670(C15orf53 Antibody) at dilution of 1:50



(under 40x lens)

A431 cells were subjected to SDS PAGE followed by western blot with Catalog No:108670(C15orf53 Antibody) at dilution of 1:1000