

# ERO1LB antibody

Catalog Number: 110370

## Product name

ERO1LB antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

ERO1LB Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, mouse brain tissue, rat brain tissue. Positive IHC detected in human pancreas tissue. Observed molecular weight by Western-blot: 54kd

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_019891). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:110370(ERO1LB antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human pancreas using Catalog No:110370(ERO1LB antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human pancreas using Catalog No:110370(ERO1LB antibody) at dilution of 1:100 (under 40x lens)