

Rat LAG3 / CD223 / Lymphocyte activation gene 3 (His Tag) recombinant protein



Catalog Number: 501370

General Information

Protein Construction

A DNA sequence encoding the rat LAG3 (Q5BK54)(Met1-Leu442) was expressed with a polyhistidine tag at the C-terminus.

Organism

Rat

Expression Host

Human Cells

QC Testing

Purity

> 90 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Gly 24

Molecular Mass

The recombinant rat LAG3 comprises 430 amino acids and predicts a molecular mass of 47.1 kDa. The apparent molecular mass of the recombinant protein is approximately 55 kDa in SDS-PAGE

under reducing conditions due to glycosylation.

Formulation

Lyophilized from sterile 20 mM NaAc, 500 mM NaCl, 10 % Glycerol, 0.02 % tween 20, pH 5.5. 1. 5 % trehalose and mannitol are added as protectants before lyophilization. 2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Rat LAG3 / CD223 / Lymphocyte activation gene 3 Protein (His Tag) SDS-PAGE