

# MAdCAM1 antibody



Catalog Number: 112370

## Product name

MAdCAM1 antibody

## Specificity

Human; other species not tested.

## Antibody description

MAdCAM1 Rabbit Polyclonal antibody. Positive IHC detected in human tonsillitis tissue, human small intestine tissue.

## Preparation

This antibody was obtained by immunization of MAdCAM1 recombinant protein (Accession Number: BC142629). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IHC

## Dilutions


Recommended Dilution:

IHC: 1:20-1:200

## Validations



Immunohistochemistry of paraffin-embedded human small intestine tissue slide using Catalog No:112370(MADCAM1 Antibody) at dilution of 1:200 (under 40x lens).



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:112370(MADCAM1 Antibody) at dilution of 1:200 (under 40x lens).