

Mouse MSR1/CD204 (His Tag) recombinant protein



Catalog Number: 502770

General Information

Gene Name Synonym

Macrophage acetylated LDL receptor I and II;
Scavenger receptor type A

Protein Construction

A DNA sequence encoding the mouse MSR1 (AAH03814.1) extracellular domain (Trp 83-Val 354) was expressed, fused with a polyhistidine tag at the N-terminus.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Purity

> 90 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

His

Molecular Mass

The recombinant mouse MSR1 consists of 288 amino acids and has a calculated molecular mass of 32 KDa.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Mouse MSR1 / CD204 / SCARA1 Protein (His Tag)