

Glutamate receptor 2 antibody



Catalog Number: 111080

Product name

Glutamate receptor 2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Glutamate receptor 2 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, human brain tissue, SKOV-3 cells. Positive IP detected in mouse brain tissue. Observed molecular weight by Western-blot: 94-105 kDa

Preparation

This antibody was obtained by immunization of Glutamate receptor 2 recombinant protein (Accession Number: NM_001083619). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:111080(GRIA2 antibody) at dilution of 1:1000



IP Result of anti-GRIA2 (IP:Catalog No:111080, 4ug; Detection:Catalog No:111080 1:500) with mouse brain tissue lysate 4000ug.