

NDUFC2 antibody



Catalog Number: 113080

Product name

NDUFC2 antibody

IHC: 1:20-1:200

Specificity

Human, Mouse; other species not tested.

Antibody description

NDUFC2 Rabbit Polyclonal antibody. Positive WB detected in human kidney tissue, A549 cells, human liver tissue, mouse brain tissue. Positive IHC detected in human kidney tissue, human hepat°Cirrhosis tissue. Observed molecular weight by Western-blot: 14kd

Preparation

This antibody was obtained by immunization of NDUFC2 recombinant protein (Accession Number: NM_004549). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



human kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:113080(NDUFC2 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human kidney using Catalog No:113080(NDUFC2 antibody) at dilution of 1:400 (under 10x lens)



Immunohistochemical of paraffin-embedded human kidney using Catalog No:113080(NDUFC2 antibody) at dilution of 1:400 (under 40x lens)