

GBE1 Polyclonal Antibody



Catalog Number: 164080

Product name

GBE1 Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to GBE1

Preparation

Antigen: A synthetic peptide of human GBE1

Formulation

PBS with 0.02% sodium azide, pH7.3.

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IF

Dilutions

WB 1:1000 - 1:2000

IF 1:20 - 1:50

Validations



Western blot - GBE1 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using GBE1 antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.