

Anti-VEGFR2/Flk-1/CD309/KDR antibody



Catalog Number: 106080

Product name

Anti-VEGFR2/Flk-1/CD309/KDR antibody

Immunogen

[Mouse VEGFR2/Flk-1/CD309/KDR \(His Tag\) recombinant protein](#)

Specificity

Mouse VEGFR2 / Flk-1 / CD309 / KDR

Antibody description

Rabbit monoclonal to VEGFR2/Flk-1/CD309/KDR

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Mouse VEGFR2 / Flk-1 / CD309 / KDR (rM VEGFR2 / Flk-1 / CD309 / KDR; P35918-1; Met1-Glu762) and conjugated with PE under optimum conditions, the unreacted PE was removed.

Formulation

Aqueous solution containing 0.5% BSA and 0.09% sodium azide

Storage

This antibody is stable for 12 months from date of receipt when stored at 2°C-8°C. Protected from prolonged exposure to light. Do not freeze ! Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

FCM

Dilutions

FCM

Validations



VEGFR2 / Flk-1 / CD309 / KDR Antibody (PE),
Rabbit MAb, Flow Cytometry

Flow cytometric analysis of Mouse KDR(CD309) expression on bEnd.3 cells. Cells were stained with PE-conjugated anti-Mouse KDR(CD309). The fluorescence histograms were derived from gated events with the forward and side light-scatter characteristics of intact cells.