

Anti-HAUSP/USP7 Rabbit antibody



Catalog Number: 167080

Product name

Rabbit IgG

Anti-HAUSP/USP7 Rabbit antibody

Applications

Specificity

WB, IHC-P, ICC/IF, FC

Human, Mouse, Rat

Dilutions

Antibody description

WB: 1/1000

HAUSP/USP7 Rabbit polyclonal antibody

IHC: 1/20

Preparation

ICC/IF: 1/50


Antigen: A synthetic peptide of human USP7

FC: 1/20

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Validations



Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Western blot detection of USP7 in C6 cell lysates using USP7 Rabbit pAb(1:1000 diluted).Predicted band size:128kDa.Observed band size:128kDa.

Ig Type