

Anti-Caspase-9 (1D1) Mouse antibody



Catalog Number: 168080

Product name

Anti-Caspase-9 (1D1) Mouse antibody

Specificity

Human

Antibody description

Caspase-9 (1D1) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human Caspase-9 protein fragments expressed in E.coli.

Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG1


Applications

WB


Dilutions

WB: 1/1000

Validations



Western blot detection of Caspase-9 in A549, MCF7, Jurkat, K562 and HeLa cell lysates using Caspase-9 mouse mAb (1:1000 diluted). Predicted band size: 49/37KDa. Observed



band size: 49KDa.

Western blot detection of Caspase-9 in CHO-K1 cell lysate (B) and CHO-K1 transfected by Caspase-9 (A) cell lysate using Caspase-9 mouse mAb (1:1000 diluted). Predicted band size: 49/37KDa. Observed band size: 49/37KDa.