

Mouse RARRES2 / TIG2 (Fc Tag) recombinant protein



Catalog Number: 503480

General Information

Gene Name Synonym

Chemerin

Protein Construction

A DNA sequence encoding the mouse Rarres2 (NP_082128.1) (Met1-Ser156) was expressed with the Fc region of human IgG1 at the C-terminus.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Purity

> 90 % as determined by SDS-PAGE.

Endotoxin

< 1.0 EU per μg protein as determined by the LAL method.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Thr 17

Molecular Mass

The recombinant mouse Rarres2 consists 378 amino acids and predicts a molecular mass of 42.6 kDa.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide

Storage


Store it under sterile conditions at -20°C to -80°C .

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Mouse RARRES2 / TIG2 Protein (Fc Tag)