

# Anti-gp130/CD130 antibody



Catalog Number: 104480

## Product name

Anti-gp130/CD130 antibody

## Immunogen

[Rat gp130/CD130 \(His Tag\) recombinant protein](#)

## Specificity

Rat gp130 / CD130 / IL6ST

## Antibody description

Rabbit monoclonal to gp130/CD130

## Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Rat gp130 / CD130 / IL6ST (rR gp130 / CD130 / IL6ST; NP\_001008725.2; Met4-Glu618) and conjugated with PE under optimum conditions, the unreacted PE was removed.

## Formulation

Aqueous solution containing 0.5% BSA and 0.09% sodium azide

## Storage

This antibody is stable for 12 months from date of receipt when stored at 2°C-8°C. Protected from prolonged exposure to light. Do not freeze ! Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG


## Applications

FCM

## Dilutions

FCM

## Validations



gp130 / CD130 / IL6ST Antibody (PE), Rabbit MAb, Flow cytometric analysis

Flow cytometric analysis of Rat IL6ST(CD130) expression on rat splenocytes. SD splenocytes were stained with PE-conjugated anti-Rat IL6ST(CD130). The histogram were derived from gated events with the forward and side light-scatter characteristics of intact cells.